

# Performance and reliability testing of a new capillary based non-contact liquid handling system

**M.Röhl (1), M. Karnath (1), U. Holfter (2), F.H. Büttner (1), R. Ries (1)**

(1) Boehringer Ingelheim Pharma GmbH & Co. KG, Birkendorfer Strasse 65, D-88397 Biberach an der Riss, Germany, www.boehringer-ingelheim.com  
 (2) CyBio AG, Göschwitzer Strasse 40, D-7745 Jena, www.cybio-ag.com

## Introduction

Direct transfer of compounds from a source plate into an assay plate can bring desirable quality improvements to screening and can increase throughput. Pre-dilution with aqueous buffers of DMSO compound solutions sometimes results in precipitation of samples which for highly active compounds may lead to false positive results if tips are not washed sufficiently. Pre-dilution is also more time consuming.

The Boehringer-Ingelheim HTS Group in Biberach Germany is looking for a reliable solution for a direct compound transfer. As a consequence of this need, an evaluation of the performance and reliability of a new capillary based non-contact liquid handling system. The CyBio®-HummingWell (Fig.1) was carried out in 2008. The CyBio®-HummingWell transfers nanolitre volumes using an array of precisely sized glass capillaries. The predefined volume to be transferred is drawn into the capillaries automatically on liquid contact and is dispensed into the source plate by compressed air pressure. The evaluation was performed in two steps. First the precision and accuracy of the CyBio®-HummingWell was tested using different cassettes (transfer volume 25 nl, 50 nl, 100 nl, 250 nl). Second, to assess stability of the device under HTS conditions, the reliability of the liquid handler was tested with 300 plates over three days under different conditions.

Figure 1: CyBio®-HummingWell device with plate stacker and wash station




Contact: Robert Ries, Boehringer Ingelheim Pharma GmbH & Co. KG  
 Dept. of Lead Discovery  
 Phone: +49 (7531) 54-4253  
 Email: robert.ries@boehringer-ingelheim.com  
This poster has been presented at MIP 150 (2008), V113-2008, Basel (October 2008)

## Material:

- Compound Storage Plate (Greiner 384 Small Volume PP Microplate, #784201)
- Assay Plate (Greiner 384 PS Transpore, #781101)
- OrangeDye, Fluka, #75380
- Dimethylsulfoxide 99.95% GC (DMSO), Sigma, #D5879-1GA
- FBS Buffer
- Multidrop Combi Dispenser, Thermo
- Flexdrop IV, Perkin Elmer
- Centrifuge, Kendro
- Saffire II Reader, Tecan
- Tecan TRAC HTS System (Figure 3)
- CyBio®-HummingWell with 2 plate stackers
- 384 Well Capillary cassettes (25nl, 50nl, 100nl, 250nl)

Figure 2: CyBio®-HummingWell capillary cassette in a 384 well plate



## Precision and accuracy test:

The testing of the precision and accuracy of the CyBio®-HummingWell was performed with 4 different capillary cassettes (25nl, 50nl, 100nl, 250nl) and with 2 different source volumes, 1µl and 2µl, in the compound plate.

The testing was done according to the following procedure:

- Filling of compound plates with 1 µl / 2 µl Orange G / DMSO solution, (Flexdrop PE)
- Visual control and centrifugation
- Aspirate (25nl, 50nl, 100nl, 250nl) Orange G / DMSO out of the compound plate (HummingWell) (automatic on immersion of capillaries)
- Dispense into the dry Assay Plate (HummingWell) (by air pressure pulse)
- Dilution of transferred Orange G / DMSO with 50µl PBS Buffer (Multidrop Combi)
- Centrifugation (1000rpm)
- Signal Readout with Saffire II (Absorbance 485nm)

Table 1: Precision an accuracy data for the 25nl, 50nl, 100nl and 250nl cassette

Cassette	25nl				50nl				100nl				250nl			
	1	2	3	4	1	2	3	4	1	2	3	4	1	2	3	4
Source Vol	2µl	1µl														
Mean (SD)	0.76 (0.40)	0.78 (0.78)	0.76 (0.71)	0.75 (0.52)	0.76 (0.48)	0.76 (0.47)	0.77 (0.47)	0.77 (0.47)	0.76 (0.48)	0.76 (0.47)	0.77 (0.47)	0.77 (0.47)	0.76 (0.48)	0.76 (0.47)	0.77 (0.47)	0.77 (0.47)
CV	0.54	0.04	0.04	0.04	0.06	0.03	0.03	0.03	0.03	0.03	0.03	0.03	0.03	0.03	0.03	0.03
Min (Max)	0.40 (0.70)	0.67 (0.69)	0.69 (0.63)	0.56 (0.41)	0.4 (0.4)	0.39 (0.39)	0.39 (0.39)	0.39 (0.39)	0.39 (0.39)	0.39 (0.39)	0.39 (0.39)	0.39 (0.39)	0.39 (0.39)	0.39 (0.39)	0.39 (0.39)	0.39 (0.39)
Max (CV)	1.97 (1.03)	1.99 (1.00)	1.92 (1.01)	1.16 (0.56)	1.97 (1.03)	1.97 (1.03)	1.97 (1.03)	1.97 (1.03)	1.97 (1.03)	1.97 (1.03)	1.97 (1.03)	1.97 (1.03)	1.97 (1.03)	1.97 (1.03)	1.97 (1.03)	1.97 (1.03)

Table 1 gives an overview of the precision and accuracy of the HummingWell device with different capillary cassettes and with 2 different compound source volumes. No outliers were observed. This shows that direct compound transfer can be performed with the HummingWell from a very small source volume of 1µl.

Figure 3: Tecan TRAC HTS System (BI Tecan 3)



## 2. Reliability of the HummingWell:

The following expanded testing routine was performed under HTS conditions, to assess the reliability of this liquid handler. The testing was done according to the following procedure:

- Compound plates (CP) filled with 2 µl Orange G / DMSO
- Aspirate Orange G / DMSO from CP (HummingWell) (automatic on immersion of capillaries)
- Dispense into the dry Assay Plate (AP) (HummingWell) (by air pressure pulse)
- Capillary Wash with 50% DMSO/Water
- Dispense 50µl Dilution buffer into AP ("Multidrop)
- Centrifugation ("VSPIN)
- Signal Readout performed Absorbance 492nm ("EnVision)
- (\* Signal Readout over night with Tecan HTS System)

Table 2: HTS Runs performed with HummingWell

Run No.	Transfer Volume	Number of MTP	Transferred Compounds
1	50nl	100	36 800
2	50nl	100	36 800
3	50nl	100	36 800
4	100nl	100	36 800
5	100nl	100	36 800
6	100nl	100	36 800

Graph 1-6: The measured signal (OD492nm) for each well is shown as a blue spot in the y-axis against the MTP number in the x-axis. Red spots are high controls, yellow spots are low controls. Pipetting error can be easily seen as blue spots below the main body of data.








Table 3: Overview of pipetting errors (outliers)

Run No.	Transfer Volume	Outliers
1	50nl	0
2	50nl	1
3	50nl	1
4	100nl	31 (26 x Capillary M14)
5	100nl	8 (4 x Capillary M14)
6	100nl	59 (57 x Capillary M14)

## Conclusion

The precision of the HummingWell is much better than the specification of 10%. The accuracy data from the 100nl cassette were higher than 10%, but still in an acceptable range. The accuracy from the 25nl and 250nl cassettes are below 10%. The robust transfer of compounds from very low source volumes (1µl) is possible with this device.

The results from the 50nl reliability test demonstrate the robust transfer of about 110.000 wells with only 2 outliers. A single capillary which often gave a lower pipetting volume was observed, but this could have been fixed easily by changing this capillary.

The data from the 100nl cassette show a different situation. One capillary (M14) did not work reproducibly. Some other capillaries failed randomly. The reasons for these outliers are not fully clear. In routine operation the M14 capillary would have been changed in this case. Finally, we believe that good quality control of the capillaries, validated cleaning and storage procedures from the vendor and intermediate function tests by the user, will all help to support the successful use of this technology in drug discovery.