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Table of Contents

1	Basic Information	17
2	Safety instructions	18
	Safety instructions for handling ortho-phosphoric acid	18
3	Technical data	19
4	Analytical method of measuring TIC	19
5	Manual TIC solids module setup	20
	5.1 TIC reaction chamber	20
	Flow setting of the carrier gas flow.....	20
	5.2 Supply of acid.....	21
	5.3 Measuring gas drying and cleaning.....	21
	5.4 Connections.....	22
6	Operation	23
	6.1 Preparing the manual TIC solids module and analyzer	23
	multi N/C 2100S and multi N/C UV HS	23
	multi N/C 3100.....	23
	6.2 Performing a TIC measurement.....	23
7	Maintenance and care	26
	7.1 Maintenance intervals	26
	7.2 Changing the desiccant in the drying tube	26
	7.3 Replacing the halogen trap	27
	7.4 Cleaning the TIC reactor	27

Index of Figures

Fig. 1	Manual TIC solids module	17
Fig. 2	TIC reactor of the manual TIC solids module	20
Fig. 3	Flowmeter for the oxygen flow.....	20
Fig. 4	Metering pump for acid on the manual TIC solids module	21
Fig. 5	Components for measuring gas drying and cleaning on the manual TIC solids module	22
Fig. 6	Connections on the manual TIC solids module.....	22

1 Basic Information

The manual TIC solids module allows the determination of TIC in solids in conjunction with the multi EA 4000 or the multi N/C 3100/2100/2100S/UVHS.

Acidification of the samples and purging of TIC is performed at elevated temperatures. The measuring gas is transported with the carrier gas to the detector of the basic unit.

The manual TIC solids module is designed exclusively for the manual supply of samples only. Automated sample supply with an autosampler is not possible.

Fig. 1 Manual TIC solids module

2 Safety instructions

IMPORTANT

This user manual is only valid in conjunction with the following documents:

- User manual of the basic module
- Software manual of the control and analysis software multiWin

Observe in particular the notes in the chapter Safety of the user manual of the basic module. The notes provided in it also apply analog without restriction to the manual TIC solids module. Special hazards emanating from the manual TIC solids module are pointed out in the chapters below.

Safety instructions for handling ortho-phosphoric acid

The ortho-phosphoric acid is used during analysis. Always observe the appropriate regulations and the notes in the EC safety data sheets of the manufacturers with regard to storage, use and disposal.

When handling ortho-phosphoric acid observe in particular the following notes:

- The concentrated acid can cause burns and irritation!
 - Always wear goggles and protective gloves when handling concentrated acids!
 - Ensure adequate ventilation in the laboratory rooms.
 - Ortho-phosphoric acid must not enter the sewage system, surface water or ground water! The applicable regulations for disposal must be meticulously observed.
 - Only store the concentrated acids in their original containers! Observe the notes on the labels.
-

IMPORTANT

In case of chemical burns from acid immediately carry out the following first aid measures:

- for skin burns:**
Immediately remove contaminated clothing and shoes, dab affected skin areas with polyethylene glycol 400 then rinse thoroughly with plenty of water.
- for eye burns:**
Rinse the eyes with open eye lid several minutes with running water (min. 10 minutes), then immediately consult a physician (even if there are no apparent symptoms).
- for burns of the respiratory organs:**
Bring the affected person into fresh air.
- if swallowed:**
Give plenty of water to drink (possibly several liters). Avoid vomiting (perforation risk). Do not attempt neutralization.

Always immediately consult a physician after carrying out the first aid measures!


3 Technical data

General characteristics	
Designation/type	Manual TIC solids module
Dimensions	300 x 500 x 550 mm
Mass	Approx. 10 kg
Procedural data	
Measuring method	TIC
Detection	NDIR (coupled with VITA method)
Sample volume	Up to 3000 mg
Gas supply	Oxygen at least 2.5 2 – 4 bar
Electrical variables	
Connection	230 V AC, optionally 115 V AC, 50/60 Hz
Protection	230 V: T10 A H, 115 V: T10 A H (Use only original fuses of Analytik Jena AG)
Typical average power consumption	415 VA
Interference suppression (electromagnetic compatibility)	in accordance with the rules of EN 55011 Group 1, Class B interference-protected (according to EN 61326-1 suitable for use in basic electromagnetic environments)

4 Analytical method of measuring TIC

During TIC analysis, the total inorganic carbon from carbonates and hydrogen carbonates in solid and powder samples is detected. Cyanides, cyanates, isocyanates and carbon particles are not detected.

To determine the inorganic carbon (TIC), acid is added to an aliquot of the sample. This leads to the decomposition of the carbonates/hydrogen carbonates contained therein. The CO₂ is purged and detected. No thermal oxidation takes place.

The measuring gas is detected in the basic device.

5 Manual TIC solids module setup

5.1 TIC reaction chamber

A 50 ml Erlenmeyer flask is used as the reactor to which a TIC head with three connections for the supply of phosphoric acid, oxygen as the carrier gas and the measuring gas outlet is attached.

A magnetic stirrer with a heating work plate ensures the sample is mixed well with the added acid and heats the mixture slightly to speed up the reaction.


Fig. 2 TIC reactor of the manual TIC solids module

Flow setting of the carrier gas flow

The carrier gas flow is set via the "O₂" needle valve with measuring glass on the back wall in the interior of the manual TIC solids module.


Fig. 3 Flowmeter for the oxygen flow

5.2 Supply of acid

The acid is metered with the metering pump by hand. The hose for supplying acid to the reactor is connected to the metering pump with a valve. The valve in the faucet of the metering pump must be closed after acid has been supplied to the reactor. This prevents acid from being inadvertently drawn into the analyzer.


Fig. 4 Metering pump for acid on the manual TIC solids module

5.3 Measuring gas drying and cleaning

A drying tube, filled with magnesium perchlorate as a desiccant, is installed in the gas path immediately downstream of the reactor.

The desiccant cartridge is followed by two water traps. The larger water trap (TC prefilter) holds back aerosols during operation, the smaller water trap (one-way retention filter) holds back rising water.

A halogen trap is used to remove interfering substances from the measuring gas and to protect the detectors and the flowmeter in the basic module. The U-tube of the halogen trap is filled with a special copper wool and brass wool. The filling of the halogen trap has to be renewed at the latest once half of the copper wool or the brass wool is discolored.


Fig. 5 Components for measuring gas drying and cleaning on the manual TIC solids module

5.4 Connections


Fig. 6 Connections on the manual TIC solids module

6 Operation

6.1 Preparing the manual TIC solids module and analyzer

The manual TIC solids module only works in conjunction with the basic module. Here, the measuring gas is not sucked through the detector with the basic unit pump, but is purged with acid in the closed system and transported with the carrier gas (oxygen) to the detector.

multi N/C 2100S and multi N/C UV HS

1. Disconnect the gas connection between the "internal" and "analyte" connections at the rear of the TOC basic unit by loosening and removing the screw connection at the "analyte" connection. (All other connections can remain as they are).
2. Connect the "analyte" connection of the manual TIC solids module to the "analyte" connection of the TOC basic unit with the enclosed tube.

multi N/C 3100

Connect the "analyte" connection of the manual TIC solids module to the "analyte" connection of the multi N/C 3100 with the enclosed tube.

CAUTION

Make sure that the "pump" connection on the TOC basic unit is not used.

6.2 Performing a TIC measurement

IMPORT

Only use dry 50 ml Erlenmeyer flasks with a standard ground joint of NS 29/32 for analysis purposes.

The sample should be pulverized as finely as possible and homogeneous.

In addition to the safety instructions for the analyzer provided in the manual of the basic unit, the following instructions also apply specifically when operating the manual solids module:

WARNING when handling phosphoric acid!

The reactor vessel and the hoses for supplying acid contain phosphoric acid! Phosphoric acid is irritating to eyes, skin and mucous membranes!

Perform slow metering operations with the dispenser to avoid overpressure!

Always wear goggles and protective gloves when handling concentrated phosphoric acid!

Rinse affected skin with water immediately.

CAUTION! Operating the magnetic stirrer!

- Do not operate the heating plate without supervision!
 - Avoid overheating the acid bath! Limit the maximum temperature to 80°C.
 - Place the glass flask in the center of the heating plate to ensure that the magnet does not hit against the wall of the flask.
 - Do not set the stirring speed of the magnet too high. The sample must not be spun against the wall of the flask.
-

1. Fill the acid bottle with 30 – 40% phosphoric acid and place it into the plastic insert in the upper housing.
2. Screw the dispenser onto the bottle. Connect the hose for supplying acid to the valve and to the TIC head.
3. Switch on both the TOC basic unit and the manual TIC solids module at the main switch.
 - ✓ **The status light on the front of the device lights up.**
4. Switch on the heating plate.

Adjust the temperature using the left rotary switch on the heating plate. A temperature setting of approx. 80°C is recommended. If the temperature selected is too high, too much water will evaporate; if the temperature is too low, the reaction will take longer.
5. Start the multiWin program.
6. Select CONFIGURATION / OPTIONS / ANALYZER COMPONENTS and activate the EXTERNAL SOLIDS MODULE option.
7. Use the METHOD / NEW menu command to create a new method.
8. Activate the HORIZONTAL FURNACE and IC measurement options.
9. Set the furnace temperature to 0°C.
10. Adjust the oxygen flow at the valve to 16 l/h.
11. Weigh the finely pulverized sample into the Erlenmeyer flask.

Note: When transferring the sample, ensure that it does not adhere to the wall of the flask. Only the bottom of the Erlenmeyer flask is covered with acid. Any sample material adhering to the walls will not be included in the analysis.
12. Carefully place the magnetic stirring rod into the Erlenmeyer flask and position the flask in the center of the heating plate.
13. Attach the TIC head to the flask. Turn the flask and the TIC head and press both lightly against each other. The ground joint connection must be gas-tight to ensure no measuring gas is lost.
14. Start the measurement by clicking on [START MEASUREMENT].
15. Enter the SAMPLE ID and, if applicable, a name in the analysis table.
16. Define the SAMPLE TYPE for the solid sample to be measured.
17. If applicable, you can enter explanations about the measurement via the [COMMENT] button.
18. With [START ►] open the MEASUREMENT window.

19. Start the measurement by clicking on the [START F2] button and follow the instructions of the control and analysis software.
20. Enter the sample quantity in [mg] in the SAMPLE window. Observe the NDIR value in the measurement window. Wait until it has passed through a maximum and then exit the SAMPLE window with [OK].
21. Wait for a message indicating operational readiness from the analyzer.
Before the measurement is started, the analyzer must be conditioned with oxygen. For this, the air must be driven out of the flask and the hose system and only oxygen must still reach the detector. Depending on the flask size, this can take up to 2 minutes. An autozero of the detector is performed at the end.
22. After the "PLEASE PRESS [OK] TO START INTEGRATION AND THEN ADD ACID!" message appears, click on [OK].
First meter 2 ml of phosphoric acid.
Slowly adjust the stirrer to an appropriate stirring speed. The sample must not be spun against the wall of the flask. Avoid the stirrer from moving in a jerky or rapid manner. Add 2 ml of acid 1 – 3 times in each case, evenly depending on the concentration and the sample quantity, until no visible reaction can be discerned anymore (foaming, formation of gas bubbles).

IMPORTANT

Since the release of CO₂ is also dependent on the metering speed and the stirring speed of the stirrer magnet, always ensure the same cycle and the same speeds are used in each measurement.

23. Close the faucet on the acid dispenser immediately after the addition of acid.
24. At the end of the measurements, the results will appear in the analysis report or in the selected analysis table.
25. At the end of the measurement, stop the stirrer magnet and open the flask:
Detach the TIC head with a twisting motion. Move the arm of the stand up, swivel it backwards and lock it in this position.
26. Clean and dry the flask.
27. Connect a dry Erlenmeyer flask to the TIC head during the measurement periods. This maintains the oxygen flush (conditioning) and the dry gas flow frees the hose system of any moisture residues.
28. Keep the heating plate at a constant temperature. Only switch off the thermostat during longer measuring periods.


7 Maintenance and care

7.1 Maintenance intervals

Maintenance task	Maintenance interval
Remove splashes of chemicals from the housing	Immediately
Clean the TIC reactor (Erlenmeyer flask)	After each sample
Change the desiccant in the drying tube	As required (blackening or clumping due to moisture absorption)
Replace the halide trap	When half of the copper wool or the brass wool is discolored
Replace the water traps	As required, but no later than after 3 months

7.2 Changing the desiccant in the drying tube

The desiccant must be changed if water absorption is visible. If the desiccant is discolored by rust or other particles, it has to be changed immediately.

1; 6 Screw closures
2; 5 Retainer clamps
3 Desiccant (magnesium perchlorate)
4 Quartz wool

1. Loosen the screw closures on the glass tube.

Caution:
The seals must not be lost and must be reinstalled in the same manner.

2. Remove the desiccant and the quartz wool.
3. Rinse the glass tube and dry it.
4. Seal the lower end of the glass tube with approx. 2 cm of quartz wool.
5. Fill the glass tube with 50 – 60 g of fresh desiccant (magnesium perchlorate).
6. Firmly screw the glass tube and push it back into the holder. Be careful not to kink the hoses.


7.3 Replacing the halogen trap

CAUTION

Risk of damage to optical and electronic components (detectors, flow sensors) from aggressive combustion products if the copper wool in the halogen trap is used up!

Replace the complete filling of the halogen trap as soon as half of the copper wool or the brass wool is discolored.

- | | |
|---|-------------|
| 1 | Gas inlet |
| 2 | Gas outlet |
| 3 | Copper wool |
| 4 | Brass wool |
| 5 | Clamp |

1. Pull the U-tube from the holder.
2. Unscrew the red union nuts from the U-tube and pull the U-tube out of the clamps (3).
3. Pull out the depleted copper wool or brass wool from the U-tube with a tweezers or a small hook.
4. Check the U-tube for cracks.

Note:
Only reuse a fully intact U-tube!
5. If required, flush the U-tube with ultrapure water and leave to dry completely.
6. Fill the U-tube with new copper wool and brass wool using tweezers or a small hook.

Note:
Replace the complete contents of the U-tube. When filling the halogen trap, make sure that the copper wool and brass wool is not compacted too much and no larger empty spaces are created in the U-tube.
7. Cover the copper wool and the brass wool with cotton wool.
8. Press the filled U-tube carefully into the clamps.
9. Close the gas inlet on the branch with copper wool and the gas output on the branch with brass wool.

7.4 Cleaning the TIC reactor

After each measurement, rinse the TIC reactor (Erlenmeyer flask) with ultrapure water. Dry it thoroughly afterwards.

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© Copyright 2017, Analytik Jena AG