# Safety manual

## PlasmaQuant 9200 (Elite)

En-us  
De  
Fr  
Es  
Bg  
Da

Et  
Fi  
El  
It  
Lv  
Lt

NI  
Pl  
Pt  
Ro  
Sv  
Sk

SI  
Cs  
Hu

---

Manufacturer                   Analytik Jena GmbH+Co. KG  
Konrad-Zuse-Straße 1  
07745 Jena / Germany  
Phone: +49 3641 77 70  
Fax: +49 3641 77 9279  
Email: info@analytik-jena.com

Technical Service               Analytik Jena GmbH+Co. KG  
Konrad-Zuse-Straße 1  
07745 Jena / Germany  
Phone: +49 3641 77 7407  
Fax: +49 3641 77 9279  
Email: service@analytik-jena.com

To ensure correct and safe use, follow these instructions. Keep these instructions for later reference.

General information           <http://www.analytik-jena.com>

Documentation number       /

Edition                         A (06/2025)

Technical documentation     Analytik Jena GmbH+Co. KG

## Table of contents

For setup sketches and setup instructions, see the end of the document [▶ 125]

US-English.....	5
Deutsch .....	11
Français.....	17
Español.....	23
Bulgarian.....	29
Dansk.....	35
Estonian .....	40
Suomi .....	45
Ελληνικά.....	50
Italiano .....	56
Latvian .....	62
Lithuanian.....	67
Nederlands.....	72
Polski.....	78
Portuguese.....	84
Român.....	90
Svensk .....	96
Slovenská .....	101
Slovenian.....	107
Česky .....	113
Magyar .....	119
Setup sketches .....	125

# 1 Basic information

## 1.1 About this manual

This safety manual describes the following device models:

- PlasmaQuant 9200
- PlasmaQuant 9200 Elite

The models are combined in the following text. Any differences between the models are explained in the relevant section.

The device is intended for operation by qualified and trained specialist personnel in accordance with the safety manual and the operating manual.

The safety manual provides information on safety instructions that must be observed when handling the device.

Detailed information about the design and operation of the device is provided in the operating manual. The operating manual provides operating personnel with the

necessary knowledge for safe handling of the device and its components. Furthermore, the operating manual includes information on the maintenance and servicing of the device as well as information on potential causes of malfunctions and their correction.

- The printed operating manual is included with the device (in German and English).
- The maintenance and servicing instructions, including videos and animations, are integrated into the software help. You can open the software help in the ASpect PQ software via the ? | **Help topics** menu item.
- The operating manual is available in German and other languages as a PDF document on the installation CD. You can also find the operating manual on the Analytik Jena homepage in the download area of the device.

## 1.2 Intended use

The optical emission spectrometer with inductively coupled plasma (ICP-OES) is used in a chemical analysis laboratory for the analysis of liquid, mainly aqueous samples, to determine concentrations of up to 75 elements up to the trace range.

The device may only be used for the processes described in the user manual. Only the specified use is regarded to be the intended use. Using the device for any other purpose may compromise the safety of the user and the device.

The device is not suitable for solutions containing hydrofluoric acid if the nebulizer or spray chamber are made of glass or quartz. Use hydrofluoric acid-resistant components for this purpose. Special provisions shall be made for operations involving organic solvents. In addition

to apparatus-related and methodical aspects, fire and health protection for the particular organic solvent must be observed here.

High-energy electric fields are used to operate the device. The device may not be operated in an explosive environment. The responsible body at the installation site must ensure that the necessary safety precautions are taken when working with flammable or explosive samples.

Consult Analytik Jena before use if you are unsure about the potential hazards of a sample solution.

Modifications, conversions and extensions are only permitted after consultation with Analytik Jena. The responsible body at the installation site alone is liable for damage caused by unauthorized modifications, conversions and extensions.

## 2 Safety

For your own safety and to ensure error-free and safe operation of the device, please read this chapter carefully before commissioning.

Observe all safety instructions listed in this manual, as well as all messages and instructions displayed by the control- and analysis software on the monitor.

### 2.1 Safety labeling on the device

Warning and mandatory action labels have been attached to the device and must always be observed. Damaged or missing warning and mandatory action labels can cause incorrect actions leading to personal injury or material damage.

- Do not remove the warning and mandatory action signs.
- Replace damaged signs.

The following warning and mandatory action signs are used:

Warning/mandatory sign	Meaning
	General warning sign
	Warning against hot surface

Warning/mandatory sign	Meaning
	Warning of hand injuries
	Observe the operating manual.
	Disconnect the power supply before opening the device cover.
	The device contains controlled substances. Analytik Jena warrants that these substances will not be released from the device within the next 25 years provided the device is employed as intended.

### 2.2 Requirements for the operating personnel

The device must only be operated by qualified specialist personnel instructed in the use of the device. The operating personnel must meet the following requirements:

- Operate the device only after instruction and training.
- Know and avoid dangers when working with the device.
- Wear personal protective equipment such as protective gloves, lab coat and safety goggles.
- Training by Analytik Jena is recommended.

The responsible body at the installation site is accountable for compliance with safety and occupational health regulations. The responsible body must meet the following requirements:

- Provide information about national regulations on work safety and accident prevention and observe them during operation of the device.
- Instruct the operating personnel in the safe operation of the device. In doing so, also convey the contents of the manuals for the device system.

### 2.3 Safety instructions for transport and commissioning

#### Transport

The device is heavy and can tip over. There is a risk of injury when lifting and carrying, especially from unsecured parts.

- Empty the device. Secure the device components according to the instructions. Remove loose parts and pack them separately. Close the plasma compartment door.
- Only transport the device in its original packaging. Insert all transport locks.

- Use a pallet truck or other suitable lifting equipment such as a crane to transport the device.
- Only lift the device with four people. Position yourself on opposite sides of the device and hold it by the four permanently screwed-in transport handles.
- Decontaminate the device before returning it to the manufacturer. Document the cleaning measures in the decontamination record. The decontamination record is provided by customer service when a return is registered.

#### Ambient conditions during commissioning

The device is dangerous if it is installed in an unsuitable environment. If the device is installed in an unsuitable environment, its service life will be reduced, e.g., due to corrosion.

- Design the installation site in accordance to the requirements in the installation conditions with installation plan.
- The device must not be installed in an explosive environment.
- Only install the device on tabletops that are suitable for the load (at least 150 kg).
- Ensure free access to the main switch on the rear of the device and to the manual plasma deactivation switch on the right-hand side of the housing.

## 2.4 Safety instructions for operation

#### Safety circuits

The plasma compartment is shielded by the housing and the UV and EMC protective glass, as well as the metal grille on the plasma compartment door in such a way that high-frequency radiation and UV radiation are reduced to a safe level. At the same time, the plasma compartment can still be accessed for maintenance purposes.

To ensure safe plasma operation, the device monitors the following conditions using safety circuits.

- The plasma compartment door is closed.

- Keep the ventilation slits clear.

#### Electrical conditions

The device may be dangerous if the conditions for the electrical connection are not met.

- Only the Analytik Jena customer service or trained and authorized specialist personnel are allowed to install and commission the device and its system components. Unauthorized assembly and installation is not permitted.
- Only use the supplied power cable or a cable of the same size with a protective earth conductor. Do not use an extension in the supply cable.
- Connect the power plug to a proper power outlet to ensure that the device meets protection class I (ground connector). Do not invalidate the protective effect by the use of an extension line which does not have a protective conductor.
- Check the electrical requirements of the device before connecting it to the mains.
- Only connect the device and its system components to the power grid when they are switched off.
- Only connect and disconnect the connection cables between the device and its system components when the device is switched off.

- The plasma torch is in working position.
- Sufficient cooling is supplied.
- Exhaust air extraction is active.
- Argon supply is ensured.

Do not ignite the plasma unless all conditions are met. If one of the safety circuits reports a fault during operation, the device extinguishes the plasma.

- Never bypass the safety circuits.
- Work safely regardless of the function of the safety circuits: Always extinguish the plasma via the control software before opening the plasma compartment door.

### Electrical hazard

Lethal voltages may occur in the device.

The induction coil operates under high voltage. An electrical discharge may occur across larger distances and cause fatal injuries, electrical shock and skin injuries.

- Before each start-up, make sure that the device and its safety devices are in proper working condition.
- In case of malfunctions of electrical components, switch off the device immediately and disconnect it from the electrical power supply.
- Do not remove or bypass any protective devices such as the housing.
- Prevent liquid from penetrating the device.
- Do not leave any movable, electrically conductive objects in the plasma compartment. This could cause a short circuit.

### Danger due to high-frequency radiation

The plasma emits electromagnetic radiation and UV radiation. The high-frequency radiation can cause serious skin and eye damage.

- Extinguish the plasma via the control software before opening the plasma compartment door.

### Thermal hazard

The plasma is very hot (up to 10,000 °C). Contact with the plasma torch shortly after operation will cause skin burns.

Flammable materials in the plasma compartment may ignite and cause a fire.

- Allow the plasma torch and surrounding area to cool for at least 5 minutes before touching.
- Do not leave any flammable materials in the plasma compartment.

The exhaust air is hot. The device's exhaust air connection is protected from contact by its plastic front and plastic panel. However, there is a risk of burns at the exhaust air connection in the laboratory.

- Do not touch the hot exhaust air connection during operation.

### Mechanical hazard

There is a risk of crushing at the hose pump during operation. Long hair and loose clothing can become caught in the rotating pump and be pulled in.

- Maintain a safe distance from moving parts.
- Wear hair protection and close-fitting clothing.

- Only install or replace pump hoses when the hose pump is at a standstill.

The device is designed so that the plasma compartment protrudes beyond the sampling compartment with the sample feed system. When working on the sample feed system, there is a risk of bumping your head on the plasma compartment or the open plasma compartment door.

- Be aware of the risk of impact when working on the sample feed system. It is best to carry out work while seated.

The water connections are under high pressure. There is a risk of injury if a hose bursts due to material fatigue. The cooling water can leak into the laboratory and cause a slipping hazard. The biocide contained in the cooling water is a hazardous substance that can cause damage to health.

- Inspect water hoses and screw connections weekly for leaks and damage.

The glass parts can break. There is a risk of injury from broken glass and glass splinters.

- Handle glass parts with care.

### Hazard from substances

The plasma leads to the formation of ozone, nitrous gases, and toxic vapors, which can cause breathing difficulties.

- Only operate the device with active extraction.

The responsible body at the installation site is in charge of monitoring air pollution emissions and minimizing harmful effects on the environment.

The device can be used to handle hazardous substances. The responsible body at the installation site is accountable for the safe handling and disposal of hazardous substances.

- Handle samples containing hydrofluoric acid with particular care.
- When working with samples dissolved in organic solvents, take the appropriate fire safety and health protection measures.
- If the device has been contaminated with hazardous substances, decontaminate it as described in the operating manual. Before using a cleaning or decontamination procedure other than that prescribed by the manufacturer, the user is required to check with Analytik Jena that the intended procedure will not damage the device.

- Handle the following operating and auxiliary materials with particular care:

Hazardous substance	Use
Methanol	Cleaning the nebulizer
Aqua regia	Cleaning the torch
Cooling water additive	Additive for cooling water to prevent corrosion and algae growth in the cooling system

- The waste container supplied is made of polyethylene and is resistant to aqueous waste solutions. When examining organic samples, check the resistance of the waste container to the solvent used. If in doubt, use a glass waste container.

### Ergonomic hazard

The device has a modular design. Maintenance and wear parts such as the water filter are easy to access so that the ergonomic strain on operating personnel is low.

- The device must be installed in such a way that allows easy access from all sides.

The overhanging plasma compartment hangs over the sampling compartment. The relatively dark working environment can cause strain on the eyes of operating personnel over long periods.

- Ensure good local lighting.
- Switch on the optional sampling compartment lighting.

### EMC compatibility

The device has been tested for interference immunity and interference emission in accordance with EN IEC 61326-1.

- The device meets the requirements for interference immunity in accordance with Table 2 (industrial environment).
- The device meets the requirements for interference emission in accordance with EN 55011 Group 1, Class A.

Electrical and electronic assemblies are shielded from interference emission by metal components.

Strong interference from other devices may disrupt the image of the plasma camera. A visual inspection of the plasma on the device is therefore recommended.

## 2.5 Safety instructions for the operation of compressed gas containers and compressed gas systems

The operating gases are taken from compressed gas containers or local compressed gas systems. The responsible body at the installation site is in charge of ensuring the safety of the gas system.

Among other things, the responsible body at the installation site must ensure that the connection type used on the outlet side of the gas pressure regulator complies with the applicable national regulations.

For safe operation, pay particular attention to the following points: Any leaks in gas systems and pipes, except for oxygen and air, can lead to an oxygen-deficient atmosphere. There is a risk of suffocation in unventilated rooms.

If oxygen accumulates in the atmosphere, flammable substances can ignite very easily and burn violently. Contact between oil or grease and oxygen under high pressure can cause explosions.

- Work on pressurized gas vessels and systems must only be carried out by trained and competent personnel. Do not carry out any unauthorized assembly or installation work.
- Only use pressure hoses and pressure reducers for the gases for which they are intended.
- Keep pipes, hoses, screw connections, and pressure reducers for oxygen free of oil and grease.
- Inspect all pipes, hoses, and screw connections regularly for leaks and damage. Repair leaks and damage without delay.
- Shut off the gas supply to the device before working on the pressurized gas vessels or the pressurized gas system. Only restart the device after completing the work and checking that it is functioning properly.

## 2.6 Safety instructions for maintenance and cleaning

There is a risk of electric shock if contact is made with live components, which may lead to serious injury.

Unauthorized servicing can lead to maladjustment or damage of the device and its system components.

- Work on electrical components inside the device may be carried out by customer service only.
- Cladding parts may only be removed by customer service, not by the customer.
- Only carry out the maintenance actions listed in the operating manual.
- Switch off the device before maintenance and cleaning. Only work on a switched-on device if this is expressly required by the operating manual.
- Before performing maintenance work inside the device, disconnect the power plug from the power outlet.

- Shut off the gas supply to the device before performing maintenance. Only leave the gas supply open if expressly required by the operating instructions.
- Use only original spare parts, wear parts and consumables. These have been tested and ensure safe operation.
- The website of Analytik Jena provides a list of consumables available for the PlasmaQuant 9200 device.
- After maintenance, ensure that all safety devices are fully functional again.
- Clean the device with a damp, non-dripping cloth. Do not use organic solvents, abrasive cleaners or bleach.

## 2.7 Behavior during emergencies

- In dangerous situations or in the event of accidents, immediately extinguish the plasma using the manual plasma deactivation switch on the right-hand side of the housing.
- If possible, only switch off the device at the power switch after allowing a cooling down period of 30 seconds. Then disconnect the mains plugs of the device and the system components from their power sockets.
- After switching off the device, close the gas supply.

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